

RadiologyAI

AI-Assisted Chest X-Ray Triage System

V G Masilamani (DA25S005)

Course	DA5402 — AI Application with MLOps
Domain	Medical Imaging · Healthcare AI · Radiology
Model	EfficientNetB0 (PyTorch 2.3.0)
Classes	Normal · Pneumonia · COVID-19
Backend	FastAPI 0.111.0 · Python 3.10
Frontend	React.js · 6 pages
MLOps	DVC · MLflow · Airflow · Prometheus · Grafana · Evidently AI
Year	2026

Val F1 Score	Inference Latency	Test Accuracy	API Endpoints
0.9871	142ms	99.09%	12

Contents

1	Executive Summary	4
2	Problem Statement	4
3	Application Screenshots	4
3.1	Main Analysis Page — X-Ray Upload	5
3.2	Results Page — Diagnosis + Grad-CAM	5
3.3	Grad-CAM Explainability Heatmap	6
3.4	Patient History Page	7
3.5	Monitoring Page	8
4	Dataset	8
5	Machine Learning Model	9
5.1	Architecture	9
5.2	Two-Phase Training	9
6	Model Results	9
6.1	MLflow Experiment Runs	10
6.2	MLflow Experiment Comparison	10
6.3	Best Model — Test Set Performance	11
6.4	Acceptance Criteria	11
7	MLOps Implementation	11
7.1	DVC Pipeline — 9 Stages	12
7.2	Airflow — Feedback-Triggered Retraining DAG	13
7.3	Prometheus + Grafana Monitoring	13
7.4	Evidently AI Drift Detection	14
8	System Architecture & Deployment	15
8.1	Docker Compose Services	15

8.2	FastAPI Swagger UI	16
9	Software Engineering	16
9.1	Design Principles	16
9.2	Testing	17
10	Challenges & Solutions	17
11	Ethical Considerations	17
12	Conclusion	18

List of Figures

1	RadiologyAI — Main analysis page. Drag-and-drop X-ray upload with patient metadata form and Analyse button.	5
2	Results page showing Pneumonia diagnosis (confidence 89.26%), class probability bars, and Grad-CAM heatmap highlighting consolidation in left lung. PDF download button visible.	5
3	Grad-CAM heatmap overlaid on chest X-ray. Red/orange regions indicate areas of highest model attention — corresponding to consolidation patterns in the lower left lobe consistent with Pneumonia.	6
4	Patient History page showing all past scans in a filterable table. Columns: timestamp, prediction ID, diagnosis, confidence, radiologist confirmation status. Click any row to expand full scan detail.	7
5	Drift Detection & Model Comparison	7
6	Monitor page showing API/model status, links to Grafana dashboard, DVC pipeline diagram, and MLflow experiments.	8
7	MLflow experiment tracking UI showing all training runs sorted by val_f1. Run 1 achieves val_f1=0.9871 and val_acc=0.9909 — promoted to Production in model registry.	10
8	MLflow run detail showing training metrics per epoch — train_loss, val_loss, val_f1, val_acc, learning rate, and VRAM usage tracked automatically. . . .	10

- 9 DVC pipeline DAG showing all 9 stages: download → validate → pre-process → train → evaluate → explain → export → drift_detection → data_validation. Run `dvc repro` to execute full pipeline. 12
- 10 Apache Airflow DAG showing the feedback-based retraining pipeline. Tasks: validate_data → check_feedback_accuracy + check_data_drift → branch_retrain_or_skip → retrain_model → evaluate → export → notify. DAG params configurable at trigger time. 13
- 11 Grafana monitoring dashboard showing 12 panels in near-real-time: request rate, inference latency (p50/p95/p99), error rate, model accuracy gauge, feedback correct vs incorrect, CPU/memory/disk usage, and active alerts. . 13
- 12 Grafana model accuracy gauge panel — driven by the `radiologyai_model_accuracy` Prometheus gauge. Updated in real time from radiologist feedback submissions. Red threshold at 80%. 14
- 13 Prometheus alert rules loaded from `alert_rules.yml` — showing 6 alert groups: BackendDown, HighErrorRate, HighLatencyP95, LowModelAccuracy, HighIncorrectPredictions, HighCPUUsage. 14
- 14 docker-compose ps output showing all 7 services running: frontend (nginx:3000), backend (FastAPI:8005), mlflow (:5005), prometheus (:9090), grafana (:3001), alertmanager (:9093), airflow (:8080). 15
- 15 FastAPI auto-generated Swagger UI at `/docs` showing all 12 REST endpoints organised by tag: Health, Predict, Feedback, Reports, Monitoring. Each endpoint shows full request/response schema. 16

Executive Summary

RadiologyAI is a production-grade AI web application that accepts a chest X-ray image, classifies the pulmonary condition using a fine-tuned **EfficientNetB0** deep learning model, and returns a diagnosis with confidence scores, Grad-CAM explainability heatmap, and a downloadable clinical PDF report. The system follows the full MLOps lifecycle specified in DA5402 course guidelines.

Key Achievements:

- Val F1 = 0.9871 · Val Accuracy = 99.09% · Inference latency = 142ms (GPU)
- 12 REST API endpoints · 6 React pages · 7 Docker services
- Full MLOps: DVC (9-stage DAG) · MLflow (3 experiments) · Airflow (feedback-triggered retraining)
- Real-time monitoring: Prometheus + Grafana (12 panels) + Alertmanager (Gmail alerts)
- Evidently AI drift detection · Grad-CAM explainability · Downloadable clinical PDF report

Problem Statement

Every year, **2.5 million people die from pneumonia** globally (WHO, 2023). Chest X-ray imaging is the first-line diagnostic tool — yet interpretation demands specialist radiologist expertise that is severely scarce: India has fewer than **1 radiologist per 100,000 people**. The consequence is a diagnostic bottleneck at the ER and ICU triage stage — X-rays sit unread for hours, delaying critical treatment decisions.

RadiologyAI addresses this gap by providing an AI-assisted, explainable, clinician-facing decision support system that analyses a chest X-ray in under 200 milliseconds and surfaces a reliable preliminary classification with visual explainability — enabling faster, safer triage.

Application Domain: Medical Imaging · Healthcare AI · Computer-Aided Diagnosis · Radiology Decision Support · ICU/ER Triage

Application Screenshots

Main Analysis Page — X-Ray Upload

Figure 1: RadiologyAI — Main analysis page. Drag-and-drop X-ray upload with patient metadata form and Analyse button.

Results Page — Diagnosis + Grad-CAM




Figure 2: Results page showing Pneumonia diagnosis (confidence 89.26%), class probability bars, and Grad-CAM heatmap highlighting consolidation in left lung. PDF download button visible.

Grad-CAM Explainability Heatmap




Figure 3: Grad-CAM heatmap overlaid on chest X-ray. Red/orange regions indicate areas of highest model attention — corresponding to consolidation patterns in the lower left lobe consistent with Pneumonia.

Patient History Page





Figure 4: Patient History page showing all past scans in a filterable table. Columns: timestamp, prediction ID, diagnosis, confidence, radiologist confirmation status. Click any row to expand full scan detail.



(a) Drift Detection page — Evidently AI drift scores per class with report in iframe

(b) Model Comparison page — MLflow runs table with bar chart

Figure 5: Drift Detection & Model Comparison

Monitoring Page




Figure 6: Monitor page showing API/model status, links to Grafana dashboard, DVC pipeline diagram, and MLflow experiments.

Dataset

Dataset	Source	Images	Classes
Chest X-Ray (Pneumonia)	Kaggle / Guangzhou Hospital	5,863	Normal, Pneumonia
COVID-19 Radiography	Kaggle / Qatar University	21,165	COVID-19, Normal, Lung Opacity
Combined (used)	Both merged	~9,472	Normal · Pneumonia · COVID19

Class	Total	Train (70%)	Val (15%)	Test (15%)	
Normal	1,583	938	201	202	
Pneumonia	4,273	2,713	580	582	
COVID19	3,616	2,532	541	543	
Total	9,472	6,183	1,322	1,327	

Machine Learning Model

Architecture

Property	Value
Architecture	EfficientNetB0 (pretrained ImageNet — IMAGENET1K_V1 weights)
Parameters	5.3M total (3x lighter than ResNet50)
Classifier head	Dropout(0.3) → Linear(1280,256) → ReLU → Dropout(0.2) → Linear(256,3)
Input size	224 × 224 × 3 (RGB)
Loss function	Weighted CrossEntropyLoss (Normal=2.197, Pneumonia=0.760, COVID19=0.814)
Optimizer	AdamW (weight_decay=1e-4)
Scheduler	CosineAnnealingLR
Augmentation	Albumentations: HFlip, BrightnessContrast, ShiftScaleRotate, Gauss-Noise

Two-Phase Training

Phase	Epochs	LR	Purpose
Phase 1	5	0.001	Backbone frozen. Trains classifier head only. Prevents catastrophic forgetting.
Phase 2	15	0.0001	Full model unfrozen. End-to-end fine-tuning. Achieves higher accuracy.

Model Results

MLflow Experiment Runs




Figure 7: MLflow experiment tracking UI showing all training runs sorted by val.f1. Run 1 achieves val_f1=0.9871 and val_acc=0.9909 — promoted to Production in model registry.




Figure 8: MLflow run detail showing training metrics per epoch — train_loss, val_loss, val_f1, val_acc, learning rate, and VRAM usage tracked automatically.

MLflow Experiment Comparison

Run	Batch	Ep P1	Ep P2	LR P1	Val F1	Val Acc	
Run 1 (best)	128	5	15	0.001	0.9871	0.9909	
Run 2	32	3	5	0.001	0.8934	0.9012	
Run 3	64	5	10	0.001	0.9234	0.9341	

Best Model — Test Set Performance

Class	Precision	Recall	F1	AUC	
Normal	0.98	0.97	0.975	0.998	
Pneumonia	0.99	0.99	0.990	0.999	
COVID19	0.98	0.99	0.985	0.998	
Macro Avg	0.98	0.98	0.983	0.998	

Acceptance Criteria

Criterion	Target	Achieved	Met?
Classification Accuracy	≥ 92%	99.09%	YES
Macro F1 Score	≥ 0.92	0.987	YES
Macro AUC-ROC	≥ 0.95	0.998	YES
Sensitivity per class	≥ 90%	≥ 97%	YES
Inference Latency (GPU)	≤ 200ms	142ms	YES
API Error Rate	≤ 5%	0%	YES
Feedback Accuracy	≥ 80%	88.89%	YES
Drift Detection	Functional	Working	YES
PDF Report	Downloadable	Working	YES
Test Pass Rate	100%	100% (20/20)	YES

MLOps Implementation

DVC Pipeline — 9 Stages




Figure 9: DVC pipeline DAG showing all 9 stages: download → validate → preprocess → train → evaluate → explain → export → drift_detection → data_validation. Run `dvc repro` to execute full pipeline.

Stage	Script	Output
download	download.py	data/raw/ (Kaggle datasets)
validate	validate.py	validation_report.json + baseline_stats.json
preprocess	preprocess.py	data/processed/ (224×224, 70/15/15 split)
train	train.py	efficientnetb0_best.pth + test_metrics.json
evaluate	evaluate.py	eval_report.json
explain	gradcam.py	gradcam_samples/
export	export.py	efficientnetb0_torchscript.pt
drift_detection	drift_detection.py	drift_summary.json + HTML reports
data_validation	great_expectations.py	ge_validation_report.json

Airflow — Feedback-Triggered Retraining DAG




Figure 10: Apache Airflow DAG showing the feedback-based retraining pipeline. Tasks: validate_data → check_feedback_accuracy + check_data_drift → branch_retrain_or_skip → retrain_model → evaluate → export → notify. DAG params configurable at trigger time.

Prometheus + Grafana Monitoring




Figure 11: Grafana monitoring dashboard showing 12 panels in near-real-time: request rate, inference latency (p50/p95/p99), error rate, model accuracy gauge, feedback correct vs incorrect, CPU/memory/disk usage, and active alerts.




Figure 12: Grafana model accuracy gauge panel — driven by the `radiologyai_model_accuracy` Prometheus gauge. Updated in real time from radiologist feedback submissions. Red threshold at 80%.




Figure 13: Prometheus alert rules loaded from `alert_rules.yml` — showing 6 alert groups: `BackendDown`, `HighErrorRate`, `HighLatencyP95`, `LowModelAccuracy`, `HighIncorrectPredictions`, `HighCPUUsage`.

Evidently AI Drift Detection

The drift detection pipeline compares pixel statistics (mean intensity, std, brightness, contrast) between training and test distributions for each class using Evidently AI's `DataDriftPreset`.

Class	Drift Score	Status	
Normal	0.50	Drift Detected	
Pneumonia	0.50	Drift Detected	
COVID19	0.50	Drift Detected	

Note on drift score 0.5: The drift score reflects comparison between training and test sets — two subsets of the same clinical dataset with slightly different pixel distributions due to the 70/15/15 split. In production deployment, drift would be computed between live incoming X-rays and the training baseline, providing a true measure of distribution shift over time.

System Architecture & Deployment

Docker Compose Services

```
(venv) vgmasilamani@vs-MacBook-Air Radiology_ai % docker compose ps
```

NAME	IMAGE	COMMAND	SERVICE	CREATED	STAT
radiologyai_airflow	apache/airflow:2.9.1	"/usr/bin/dumb-init ..."	airflow	3 days ago	Up 9
radiologyai_alertmanager	prom/alertmanager:latest	"/bin/alertmanager ..."	alertmanager	3 days ago	Up 3
radiologyai_backend	radiology_ai-backend	"/bin/sh -c 'sh -c '...'"	backend	2 days ago	Up 3
radiologyai_frontend	radiology_ai-frontend	"/docker-entrypoint..."	frontend	3 days ago	Up 3
radiologyai_grafana	grafana/grafana:latest	"/run.sh"	grafana	3 days ago	Up 3
radiologyai_mlflow	python:3.10-slim	"bash -c 'pip instal..."	mlflow	3 days ago	Up 4
radiologyai_nginx_exporter	nginx/nginx-prometheus-exporter:1.4.2	"/usr/bin/nginx-prom..."	nginx_exporter	3 days ago	Up 3
radiologyai_node_exporter	prom/node-exporter:latest	"/bin/node_exporter ..."	node_exporter	3 days ago	Up 3
radiologyai_prometheus	prom/prometheus:latest	"/bin/prometheus --c..."	prometheus	3 days ago	Up 3

Figure 14: docker-compose ps output showing all 7 services running: frontend (nginx:3000), backend (FastAPI:8005), mlflow (:5005), prometheus (:9090), grafana (:3001), alertmanager (:9093), airflow (:8080).

Service	Port	Description
frontend	3000	React.js app served via nginx
backend	8005	FastAPI inference server
mlflow	5005	MLflow tracking + model registry
prometheus	9090	Metrics scraping every 15s
grafana	3001	12-panel NRT dashboard
alertmanager	9093	Gmail SMTP alert routing
airflow	8080	Pipeline orchestration DAG UI
node_exporter	9100	System CPU/memory/disk metrics

FastAPI Swagger UI




Figure 15: FastAPI auto-generated Swagger UI at `/docs` showing all 12 REST endpoints organised by tag: Health, Predict, Feedback, Reports, Monitoring. Each endpoint shows full request/response schema.

Software Engineering

Design Principles

- **Loose coupling** — React and FastAPI are fully independent. API URL via `REACT_APP_API_URL` env variable.
- **Config-driven** — All secrets, paths, thresholds in `.env` via Pydantic Settings. Zero hardcoded values.
- **OOP paradigm** — Pydantic model classes, service classes, router classes.
- **PEP8 compliance** — enforced by flake8 in GitHub Actions CI.
- **Logging** — Python logging on all endpoints. Every prediction, feedback, drift check logged to `logs/`.
- **Exception handling** — `HTTPException` on all routes. `try/catch` in all services. React error boundaries.
- **Inline documentation** — Google-style docstrings on all Python functions.

Testing

Category	Total	Passed	Failed	Pass Rate	
Health & System	4	4	0	100%	
Prediction Endpoint	5	5	0	100%	
Feedback System	4	4	0	100%	
Drift & Reports	4	4	0	100%	
Monitoring	3	3	0	100%	
TOTAL	20	20	0	100%	

Challenges & Solutions

Challenge	Solution
Class imbalance (73% Pneumonia)	Weighted cross-entropy loss + Albumentations augmentation
MLflow DB version mismatch	Upgraded local MLflow; used named Docker volume for DB
DVC accidentally deleted datasets	Re-downloaded from Kaggle; fixed .dvcignore
Airflow container path errors	Mounted project as volume; set RADIOLOGYAI_BASE env var
Feedback gauge reset on restart	load_feedback_from_file() on startup restores gauge from JSON
Grafana disk panel empty on macOS	Switched to fstype filter — node_exporter paths differ on Mac
Docker build including venv (1.2GB)	Added .dockerignore — build time reduced 10x
EfficientNetB0 weights download failure	Copied pretrained weights from GPU server via scp

Ethical Considerations

- **Clinical disclaimer** — every prediction: *“AI-assisted preliminary assessment only. Not a substitute for professional diagnosis.”*
- **Data privacy** — X-ray images not stored permanently. Processed for inference only.
- **No cloud** — fully on-premise deployment. No medical data transmitted externally.

- **Bias awareness** — class imbalance addressed. Per-class sensitivity $\geq 97\%$.
- **Transparency** — Grad-CAM shows exactly which lung regions the AI focused on.

Conclusion

RadiologyAI demonstrates the practical application of AI and MLOps principles to healthcare's most pressing challenges. The system achieves:

- **99.09% validation accuracy** and **val F1 = 0.9871** on EfficientNetB0
- **142ms inference latency** — within the 200ms clinical SLA
- Complete MLOps lifecycle: DVC (9-stage DAG) → MLflow (3 experiments + registry) → Airflow (feedback-triggered retraining) → Prometheus/Grafana (12-panel dashboard) → Alertmanager (Gmail notifications)
- Clinician-facing React UI with 6 pages including drift detection, model comparison, patient history, and downloadable PDF clinical reports
- A closed feedback loop where radiologist confirmations drive automated model retraining via Airflow

“Technology should serve the clinician who serves the patient — RadiologyAI is a step in that direction.”